









19. Screw Neck Vials for Storage Purposes NEW

- Screw Neck Vials for storage purposes made out of 1st hydrol. class glass
- Vials with different volumes are available
- Clear and amber Vials
- Seals with different septa materials
- Barcode labelling upon request



19.1 Screw Neck Vials for Storage Purposes NEW

 Scale 1:2, full-scale s. p. 96	 Scale 1:2, full-scale s. p. 96	 Scale 1:2, full-scale s. p. 95	 Scale 1:2, full-scale s. p. 95	 Scale 1:2, full-scale s. p. 93	 Scale 1:2, full-scale s. p. 93
11 09 0210	11 09 0259	13 09 0222	13 09 0280	15 09 1703 15 09 1774	15 09 1657 15 09 1800
1.5ml Screw Neck Vial 8-425, 32 x 11.6mm, clear glass, 1 st hydrol. class	1.5ml Screw Neck Vial 8-425, 32 x 11.6mm, amber glass, 1 st hydrol. class	4ml Screw Neck Vial 13-425, 45 x 14.7mm, clear glass, 1 st hydrol. class	4ml Screw Neck Vial 13-425, 45 x 14.7mm, amber glass, 1 st hydrol. class	8ml Screw Neck Vial 15-425, 61 x 16.6mm, clear/amber glass, 1 st hydrol. class	12ml Screw Neck Vial 15-425, 66 x 18.5mm, clear/amber glass, 1 st hydrol. class
100 pcs. per PP-Box	100 pcs. per PP-Box	100 pcs. per PP-Box	100 pcs. per PP-Box	100 pcs. per carton	100 pcs. per carton

 Scale 1:2, full-scale s. p. 93	 Scale 1:2, full-scale s. p. 93	 Scale 1:2, full-scale s. p. 93	 Scale 1:2, full-scale s. p. 93	 Scale 1:2, full-scale s. p. 93	 Scale 1:2, full-scale s. p. 93
18 09 1704	20 09 1705	24 09 0589 24 09 0927	24 09 0839 24 09 0923	24 09 0402 24 09 0928	24 09 1089 24 09 1090
16ml Screw Neck Vial 18-400, 71 x 20.6mm, clear glass, 1 st hydrol. class	20ml Screw Neck Vial 20-400, 86 x 22.7mm, clear glass, 1 st hydrol. class	20ml Screw Neck Vial 24-400, 57 x 27.5mm, clear/amber glass, 1 st hydrol. class	30ml Screw Neck Vial 24-400, 72.5 x 27.5mm, clear/amber glass, 1 st hydrol. class	40ml Screw Neck Vial 24-400, 95 x 27.5mm, clear/amber glass, 1 st hydrol. class	60ml Screw Neck Vial 24-400, 140 x 27.5mm, clear/amber glass, 1 st hydrol. class
100 pcs. per carton	100 pcs. per carton	100 pcs. per PP-Box	100 pcs. per PP-Box	100 pcs. per PP-Box	100 pcs. per PP-Box

19.2 PP Screw Seals for Storage Vials

NEW



- Seals with different septa material are available
- Corresponding centre hole versions are partially available.
- Packed with 100 pieces in tamper-proof evident zip-lock PE-bags

19.2.1 PP Screw Seals ND8 for 11 09 0210 and 11 09 0259

Art. No.	Description Cap	Septa Material	Durometer	Thickn.
08 15 0654	Polypropylene Screw Cap black, closed top	Nat. Rubber red-orange/TEF transparent	60° shore A	1.3mm
08 15 1653	Polypropylene Screw Cap black, closed top	Butyl red/PTFE grey	55° shore A	1.3mm
08 15 1040	Polypropylene Screw Cap black, closed top	Silicone white/PTFE red	45° shore A	1.3mm

Packaging Unit: 100 pcs. per PE-Bag

19.2.2 PP Screw Seals ND13 for 13 09 0222 and 13 09 0280

Art. No.	Description Cap	Septa Material	Durometer	Thickn.
13 15 0439	Polypropylene Screw Cap black, closed top	Nat. Rubber red-orange/TEF transparent	60° shore A	1.3mm
13 15 1638	Polypropylene Screw Cap black, closed top	Butyl red/PTFE grey	55° shore A	1.3mm
13 15 0648	Polypropylene Screw Cap black, closed top	Silicone cream/PTFE red	55° shore A	1.5mm

Packaging Unit: 100 pcs. per PE-Bag

19.2.3 PP Screw Seals ND15 for 15 09 1703, 15 09 1774, 15 09 1657, 15 09 1800

Art. No.	Description Cap	Septa Material	Durometer	Thickn.
15 15 0793	Polypropylene Screw Cap black, closed top	Nat. Rubber red-orange/TEF transparent	60° shore A	1.3mm
15 15 1748	Polypropylene Screw Cap black, closed top	Butyl red/PTFE grey	55° shore A	1.6mm
15 15 1083	Polypropylene Screw Cap black, closed top	Silicone white/PTFE red	45° shore A	1.3mm
15 15 1938	Polypropylene Screw Cap black, 9mm centre hole	Nat. Rubber red-orange/TEF transparent	60° shore A	1.3mm
15 15 1932	Polypropylene Screw Cap black, 9mm centre hole	Butyl red/PTFE grey	55° shore A	1.6mm
15 15 1989	Polypropylene Screw Cap black, 9mm centre hole	Silicone white/PTFE red	45° shore A	1.3mm

Packaging Unit: 100 pcs. per PE-Bag

NEW
NEW
NEW

19.2.4 PP Screw Seals ND18 for 18 09 1704

Art. No.	Description Cap	Septa Material	Durometer	Thickn.
18 15 1387	Polypropylene Screw Cap black, closed top	Butyl red/PTFE grey	55° shore A	1.6mm
18 15 1132	Polypropylene Screw Cap black, closed top	Silicone blue transparent/PTFE white	45° shore A	1.7mm
18 15 2069	Polypropylene Screw Cap black, closed top	Silicone white/PTFE red	55° shore A	1.5mm

Packaging Unit: 100 pcs. per PE-Bag

19.2.5 PP Screw Seals ND20 for 20 09 1705

Art. No.	Description Cap	Septa Material	Durometer	Thickn.
20 15 1803	Polypropylene Screw Cap white, closed top	Nat. Rubber red-orange/TEF transparent	60° shore A	1.3mm
20 15 1805	Polypropylene Screw Cap white, closed top	Butyl red/PTFE grey	55° shore A	1.3mm
20 15 1804	Polypropylene Screw Cap white, closed top	Silicone white/PTFE red	45° shore A	1.3mm

Packaging Unit: 100 pcs. per PE-Bag

19.2.6 PP Screw Seals ND24 for 24 09 0589, 24 09 0927, 24 09 0839, 24 09 0923, 24 09 0402, 24 09 0928, 24 09 1089, 24 09 1090

Art. No.	Description Cap	Septa Material	Durometer	Thickn.
24 15 1395	Polypropylene Screw Cap white, closed top	Butyl red/PTFE grey	55° shore A	2.5mm
24 15 1540	Polypropylene Screw Cap white, closed top	Silicone white/PTFE beige	45° shore A	3.2mm
24 15 1007	Polypropylene Screw Cap white, closed top	PTFE/EPDM/PTFE	65° shore A	2.0mm
24 04 0841	UltraBond Seal white, closed top	Silicone natural/PTFE beige	45° shore A	3.2mm

Packaging Unit: 100 pcs. per PE-Bag